

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3-Amino-2-chloro-6-methoxypyridine - 97%

Product Number:

681679

CAS Number:


34392-85-3

Formula:

C₆H₇ClN₂O

Formula Weight:

158.59 g/mol

**TEST****Specification**

Appearance (Color)

White to Dark Brown

Appearance (Form)

Powder or Crystal or Chunk(s)

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 96.5 %

Titration by AgNO₃

96.5 - 103.5 %

(after Oxygen Combustion)

Specification: PRD.0.ZQ5.10000024262

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.