

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com Email USA: techsery@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Blue Value

Specification: PRD.1.ZQ5.10000026363

Amylose from potato - used as amylase substrate

Product Number: A0512 CAS Number: 9005-82-7 MDL: MFCD00134652 Formula: (C6H10O5)n

TEST Specification

Appearance (Color) White to Off White

Pow der Appearance (Form) Solubility (Color) Colorless

Solubility (Turbidity) Clear to Slightly Hazy

1 mg/mL, 0.05 M NaOH

Solvent Content (by GC) < 10 %

as Ethanol

Solvent Content (by GC) < 2 %

as Butanol ≥ 0.5

Recommended Retest Period

6 years

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.