

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:


Methyltriphenylphosphonium iodide, polymer-bound – 200–400 mesh, extent of labeling: ~2.0 mmol/g loading, 2 % cross-linked with divinylbenzene

Product Number:

572942

MDL:

MFCD03703043



### TEST

### Specification

Appearance (Color)

Yellow to Brown

Appearance (Form)

Powder

Loss on Drying

4 Hours at 50 Degrees Celsius

Swelling Volume mL/g (THF)

Loading (mmol/g)

Approximately 2.0 mmol/g

Specification: PRD.1.ZQ5.10000032737

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.