

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Human IgG (γ-chain specific)-R-Phycoerythrin antibody produced in goat - affinity isolated antibody, buffered aqueous solution

Product Number: **P9170**

MDL: MFCD00162448

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Turbidity)	Clear
Appearance (Color)	Pink
Appearance (Form)	Solution
Immunofluorescence (Direct): Minimum 1:16 Cells: Human Peripheral blood Lymphocytes	Conforms
Protein Content Prior to the addition of BSA 0.5 mg/mL	
A567/A280 Ratio Prior to the addition of BSA	1.5 - 5.0

Specification: PRD.2.ZQ5.1000007119

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.