

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Bionic™ Buffer, 10X Concentrate – for Northern and Southern blotting, solution

Product Number:

B6185

TEST	Specification
------	---------------

Appearance (Turbidity)

Clear

Appearance (Color)

Colorless

Appearance (Form)

Liquid

pH

8.0 - 8.5

Neat, at Room Temperature

DNase, Exonuclease Detection

None Detected

NICKase, Endonuclease Detection

None Detected

Suitability

Suitable

Suitable for use in the electrophoretic analysis of DNA.

Specification: PRD.0.ZQ5.10000083619

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.