

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification


Product Name:
Ethoxynonafluorobutane - 99.0%, mixture of two isomers

Product Number: **592781**

MDL: MFCD01862062

Formula: C₆H₅F₉O

Formula Weight: 264.1 AMU



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 99.0 %
As a Mixture of Isomers	
Water (by Karl Fischer)	≤ 100.0 ppm
Residue on Evaporation	≤ 0.0003 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.