


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Aniline-¹⁵N – 98 atom % ¹⁵N

Product Number: 488100
CAS Number: 7022-92-6
MDL: MFCD00084195
Formula: C₆H₇¹⁵N
Formula Weight: 94.12 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Clarity)	Clear
Appearance (Color)	Colorless to Slightly Yellow
Appearance (Form)	Liquid
Purity (GC)	≥ 99.0 %
Proton NMR Spectrum	Conforms
¹³ C NMR Spectrum	Conforms
¹⁵ N Enrichment	≥ 98.0 %

Specification: PRD.0.ZQ5.10000064238

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.