

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Hydroxypropyl cellulose - average M_w ~100,000, powder, 20 mesh particle size (99% through)

Product Number:


191884

CAS Number:

9004-64-2

MDL:

MFCD00132688

**TEST****Specification**

Appearance (Color)

White

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Viscosity

65 - 175 cps

C = 5%, H₂O at 25 Degrees Celsius

Specification: PRD.2.ZQ5.10000030014

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.