

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

(R)-(+)-3-Hydroxypiperidine hydrochloride 98%

Product Number:

410462

CAS Number:

198976-43-1

MDL:


MFCD00192221

Formula:

C₅H₁₁NO · HCl

Formula Weight:

137.61 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Conforms to Requirements
Solid, Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃	97.5- 102.5 %
Optical Rotation	7.0- 9.0
Concentration = 2%, Methanol	

Remarks:

Specification Date : 06/05/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.