


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Trimethyl(2,3,4,5-tetramethyl-2,4-cyclopentadien-1-yl)silane - 97%

Product Number: 463418
CAS Number: 134695-74-2
MDL: MFCD00674076
Formula: C₁₂H₂₂Si
Formula Weight: 194.39 g/mol



TEST

Specification

Appearance (Color)	Faint Yellow to Yellow and Faint Green-Yellow to Green-Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Refractive index at 20 ° C	1.484 - 1.488
Purity (GC)	≥ 96.5 %

Specification: PRD.1.ZQ5.10000007779

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.