

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Calcium - pieces, <1 cm, 99%

Product Number: 327387 Ca
CAS Number: 7440-70-2
MDL: MFCD00085314
Formula: Ca
Formula Weight: 40.08 g/mol

TEST	Specification
Appearance (Color)	Conforms to Requirements
Grey	
Appearance (Form)	Conforms to Requirements
Pieces or Pellets	
Purity (Titration by NaOH)	98.5 - 101.5 %
Size	Conforms
Typically < 1 cm	

Specification: PRD.0.ZQ5.10000034718

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.