

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)


Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Halocarbon oil 700

**Product Number:** H8898  
**CAS Number:** 9002-83-9  
**MDL:** MFCD00147324



TEST	Specification
Appearance (Clarity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Viscosity	630 - 770 cps
At 100 Degrees Fahrenheit	
Miscellaneous Assay	≤ 10 ppm
Hydrocarbon	
Miscellaneous Assay	≥ 15 min
Saturation Test	
Acidity	≤ 0.0002 mg/g
mg KOH/g	

Specification: PRD.0.ZQ5.10000012665

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.