

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Bis(hexamethylene)triamine - high purity

Product Number:**421960**

CAS Number:

143-23-7

MDL:


MFCD00010448

Formula:

C₁₂H₂₉N₃

Formula Weight:

215.38 g/mol

**TEST****Specification**

Appearance (Color)

White to Off White

Appearance (Form)

Conforms to Requirements

Chunks and/or Solid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 94.0 %

Specification: PRD.1.ZQ5.10000011146

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.