


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Trimethylolpropane propoxylate triacrylate - average $M_n \sim 644$

Product Number: 407577
CAS Number: 53879-54-2
MDL: MFCD00197950
Formula: C₃₃H₅₆O₁₂
Formula Weight: 644.79 g/mol



TEST

Specification

Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Viscosity	70 - 170 cps
At 25 Degrees Celsius	

Specification: PRD.1.ZQ5.10000033146

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.