

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Aluminum titanate – nanopowder, &lt;25 nm particle size (BET), 98.5% trace metals basis

Product Number:

634131

CAS Number:

12004-39-6

MDL:


MFCD00064749

Formula:

Al<sub>2</sub>O<sub>5</sub>Ti

Formula Weight:

181.83 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
ICP Major Analysis Confirms Al and Ti Components	Confirmed
Trace Metal Analysis >= 98.5% by XRF	Conforms
Particle Size < 25 nm (BET)	Confirmed

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.