

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

1-O-Octadecyl-2-O-methyl-sn-glycero-3-phosphorylcholine - ≥99% (TLC), waxy solid

Product Number:

O9262



MDL:

MFCD00063527

Formula:

C₂₇H₅₈NO₆P

Formula Weight:

523.73 g/mol

Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Color)

White

Appearance (Form)

Powder or Solid

Solubility

Conforms

Soluble in EtOH at 15 mg/mL

Proton NMR Spectrum

Conforms to Structure

Purity (TLC)

≥ 99 %

Specification: PRD.1.ZQ5.10000058630

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.