


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Dacarbazine - antineoplastic purine analog

Product Number: D2390
CAS Number: 4342-03-4
MDL: MFCD00057167
Formula: C₆H₁₀N₆O
Formula Weight: 182.18 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)	Off-White to Light Yellow
Appearance (Form)	Powder
Solubility (Color)	Colorless to Yellow, Pink or Dark Orange
Solubility (Turbidity)	Clear to Very Slightly Hazy
50 mg/ml, 1 M HCl	
Infrared spectrum	Conforms to Structure
Purity (TLC)	≥ 98 %

Specification: PRD.0.ZQ5.10000028895

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.