


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
BQ-123 - ≥99%, lyophilized powder

Product Number: B150
CAS Number: 136553-81-6
MDL: MFCD00133204
Formula: C₃₁H₄₂N₆O₇
Formula Weight: 610.70 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
0.25 mg/mL in H ₂ O	
Amino Acid Analysis	≥ 65 %
Peptide Content	
Purity (HPLC)	≥ 99 %

Specification: PRD.0.ZQ5.10000044490

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.