

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Trimethyl(methylcyclopentadienyl)platinum(IV) - packaged for use in deposition systems

Product Number:

697540

CAS Number:


94442-22-5

Formula:

C₉H₁₆Pt

Formula Weight:

319.30 g/mol



TEST

Specification

Appearance (Color)

White to Light Yellow to Light Orange

Appearance (Form)

Powder or Crystals or Solid or Chunks

Proton & Carbon-13 NMR Spectra

Conforms to Structure

Carbon

33.0 - 34.7 %

ICP Major Analysis

Confirmed

Confirms Platinum Component

Specification: PRD.3.ZQ5.10000051557

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.