

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Decanoic acid-1-¹³C - 99 atom % ¹³C

Product Number: 488658
CAS Number: 84600-66-8
MDL: MFCD00083932
Formula: 13CC9H20O2
Formula Weight: 173.26 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Solid
Water (by Karl Fischer)	≤ 0.50 %
Purity (GC)	≥ 99.0 %
1H NMR Spectrum	Conforms to Structure
13C NMR Spectrum	Conforms to Structure
13C Enrichment	≥ 99.0 atom %

Specification: PRD.0.ZQ5.10000102838

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

