

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Dichloro(dicyclopentadienyl)platinum(II) - 97%

Product Number: **701203**
CAS Number: 12083-92-0

Formula: C₁₀H₁₂Cl₂Pt
Formula Weight: 398.19 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Solid
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃	≥ 97.0 %
Carbon	29.1 - 31.2 %
Recommended Retest Period	-----
2 Years	

Specification: PRD.1.ZQ5.10000011159

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.