

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Astrazon Orange G

Product Number: 274216
 CAS Number: 3056-93-7
 MDL: MFCD00011947
 Formula: C₂₂H₂₃ClN₂
 Formula Weight: 350.88 g/mol



TEST	Specification
Appearance (Color)	Red
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Wavelength c = 0.01 g/l; Methanol	268 - 274 nm
Extinction Coefficient	≥ 3000
Wavelength	282 - 288 nm
Extinction Coefficient	≥ 3000
Wavelength	487 - 493 nm
Extinction Coefficient	≥ 19000
Carbon	30.0 - 60.0 %
Solubility (Turbidity) c = 10mg/ml, MeOH:H ₂ O (1:1)	Clear to Hazy
Solubility (Color)	Red-Orange

Specification: PRD.0.ZQ5.10000066824

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.