

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

5,6-Epoxy-5,6-dihydro-[1,10]phenanthroline 98%

Product Number:

578789

CAS Number:

65115-91-5

MDL:


MFCD04039972

Formula:

C₁₂H₈N₂O

Formula Weight:

196.2 g/mol



TEST	Specification
Appearance (Color)	White to Beige
Appearance (Form) Powder or Crystals	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Carbon	71.6- 75.3 %
Nitrogen	13.9- 14.6 %
Purity (HPLC)	≥97.5 %

Remarks:

Specification Date : 06/21/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.