

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
[1,1'-Bis(di-tert-butylphosphino)ferrocene]dichloropalladium(II) - 98%

Product Number: **701602**
CAS Number: 95408-45-0

Formula: C₂₆H₄₄Cl₂FeP₂Pd
Formula Weight: 651.75 g/mol
Storage Temperature: -20 °C



TEST	Specification
------	---------------

Appearance (Color)	Beige to Very Dark Beige and Brown to Very Dark Brown and Red-Brown and Very Dark Red-Brown
Appearance (Form)	Powder
Proton NMR spectra	Conforms to Structure
Phosphorus NMR Spectra	Conforms to Structure
Chlorine	10.13 - 11.63 %
Carbon	46.7 - 49.1 %
Purity	97.5 - 102.5 %
based on Carbon content	

Specification: PRD.1.ZQ5.10000021631

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.