

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Iron(III) fluoride trihydrate

Product Number: 334499
CAS Number: 15469-38-2
MDL: MFCD00149718
Formula: F₃Fe · 3H₂O
Formula Weight: 166.89 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
White to Beige to Tan to Brown to Red	
Appearance (Form)	Powder
Titration by Na ₂ S ₂ O ₃	32.3 - 34.6 %
% Iron	
ICP Major Analysis	Confirmed
Confirms Iron Component	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.