


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Methyl indole-3-carboxylate - 99%

Product Number: 395307
CAS Number: 942-24-5
MDL: MFCD00189407
Formula: C₁₀H₉NO₂
Formula Weight: 175.18 g/mol



TEST	Specification
Appearance (Color)	White to Beige to Pink
Appearance (Form)	Powder
Infrared Spectrum	Conforms to Structure
Carbon	67.5 - 69.6 %
Nitrogen	7.7 - 8.3 %
Purity (GC)	≥ 98.5 %

Specification: PRD.1.ZQ5.10000019152

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.