


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tetramethylammonium hydroxide pentahydrate - ≥97%

Product Number: T7505
CAS Number: 10424-65-4
MDL: MFCD00149566
Formula: C₄H₁₃NO·5H₂O
Formula Weight: 181.23 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Titration with HCl	97.0 - 103.0 %
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
100 mg/ml, H ₂ O	

Specification: PRD.1.ZQ5.10000040878

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.