


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3,5-Di-tert-butylcatechol – 98%

Product Number: D45800
CAS Number: 1020-31-1
MDL: MFCD00008819
Formula: C₁₄H₂₂O₂
Formula Weight: 222.32 g/mol



TEST

Specification

Appearance (Color)

White to Light Yellow and Faint Beige to
Light Beige and Faint Brown to Light Brown

Appearance (Form)

Powder or Crystals or Chunk(s) or
Powder with Chunk(s) or Crystals with Chunk(s)

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 97.5 %

Specification: PRD.1.ZQ5.10000015647

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.