

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Cyclopropyltriphenylphosphonium bromide - 98%

Product Number:

157317

CAS Number:

14114-05-7

MDL:


MFCD00011872

Formula:

C₂₁H₂₀BrP

Formula Weight:

383.26 g/mol

**TEST****Specification**

Appearance (Color)

White to Off-White

Appearance (Form)

Conforms to Requirements

Powder, Crystalline Powder, Crystals and/or Chunks

Proton NMR spectra

Conforms to Structure

Titration with AgNO₃

97.5 - 102.5 %

Specification: PRD.0.ZQ5.10000042297

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.