


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tetramethylammonium nitrate - 95%

Product Number: 425257
CAS Number: 1941-24-8
MDL: MFCD00043103
Formula: C₄H₁₂N₂O₃
Formula Weight: 136.15 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder or Crystal or Chunk(s)
Infrared spectrum	Conforms to Structure
Titration with HClO ₄	94.5 - 105.5 %
Purity (TLC)	≥ 94.5 %
Solubility (Turbidity) 25mg/ml, Methanol	Clear
Solubility (Color)	Colorless

Specification: PRD.2.ZQ5.10000005407

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.