


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Chitosan - high molecular weight

Product Number: **419419**
CAS Number: 9012-76-4
MDL: MFCD00161512



TEST	Specification
------	---------------

Appearance (Color)	Off White to Beige to Orange
Appearance (Form)	Powder or Flakes
Assay	≥ 75 %
Decetylation	
Viscosity	800 - 2000 cps
c = 1%; 1% Acetic Acid	

Specification: PRD.4.ZQ5.10000004416

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.