

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Polycaprolactone diol – average $M_n \sim 2,000$

Product Number: **189421**
CAS Number: 36890-68-3



Formula: $(C_6H_{10}O_2)_n$

| TEST | Specification |
|---------------------------------|--------------------------|
| Appearance (Color) | White |
| Appearance (Form) Waxy Solid | Conforms to Requirements |
| Infrared spectrum | Conforms to Structure |
| Hydroxyl Value mg KOH/g | 51 - 63 |
| Molecular Number by GPC | 1500 - 2500 |
| Water (by Karl Fischer) | $\leq 1.0 \%$ |

Specification: PRD.0.ZQ5.10000044818

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.