


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Polyethylene - High density, melt index 2.2 g/10 min (190 °C/2.16kg)

Product Number: 547999
CAS Number: 9002-88-4
MDL: MFCD00084423
Formula: C₂H₄



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Pellets
Infrared spectrum	Conforms to Structure
Melt Index	2.0 - 2.4 gm 10

Specification: PRD.0.ZQ5.10000103171

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

