

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-AGO2 antibody produced in rat - ~1.5 mg/mL, clone 11A9, purified immunoglobulin

Product Number:

SAB4200085

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	1.5 - 2.0 mg/ml
WB-Cell Line/Tissue Extract	Hela Cells
WB Working Conc	0.5 - 1.0 µg/ml
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000015001

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.