

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Calcium nitrate tetrahydrate – cell culture tested

Product Number:

C7652

CAS Number:

13477-34-4

MDL:


MFCD00149604

Formula:

CaN2O6 · 4H2O

Formula Weight:

236.15 g/mol



TEST

Specification

Appearance (Color)

White

Appearance (Form)

Powder or Crystals

Solubility (Turbidity)

Clear

Solubility (Color)

Colorless

at 20g plus 100 mL of Water

Complexometric EDTA

≥ 99.0 %

Cell Culture Test

Pass

Recommended Retest Period

3 Years

Specification: PRD.0.ZQ5.10000012518

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.