

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)


Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
1,10-Phenanthroline monohydrate - ACS reagent, 99%

**Product Number:** 320056  
**CAS Number:** 5144-89-8  
**MDL:** MFCD00011678  
**Formula:** C<sub>12</sub>H<sub>8</sub>N<sub>2</sub> · H<sub>2</sub>O  
**Formula Weight:** 198.22 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form) Crystals or Crystalline Powder	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Carbon	71.6 - 73.8 %
Nitrogen	13.9 - 14.3 %
Suitability	Suitable
Suitable for Use as a Redox Indicator	
Suitability	Suitable
Suitable for Use in the Determination of Iron	
Meets ACS Requirements	Current ACS Specification

Specification: PRD.1.ZQ5.10000003568

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.