


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Caffeic acid - ≥98.0% (HPLC)

Product Number: C0625
CAS Number: 331-39-5
MDL: MFCD00004392
Formula: C₉H₈O₄
Formula Weight: 180.16 g/mol



TEST

Specification

Appearance (Color)	Faint Yellow to Yellow to Brown
Appearance (Form)	Powder
Solubility (Color)	Light Yellow to Dark Yellow to Dark Brown
Solubility (Turbidity)	Clear to Very Slightly Hazy
50 mg/mL, EtOH	
¹ H NMR Spectrum	Conforms to Structure
Purity (HPLC)	≥ 98.0 %

Specification: PRD.2.ZQ5.10000009781

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.