

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

cis-Diammineplatinum(II) dichloride ≥99.9% trace metals basis

Product Number:**479306**Pt(NH₃)₂Cl₂

CAS Number:

15663-27-1

MDL:

MFCD00011623

Formula:

H₆Cl₂N₂Pt

Formula Weight:

300.05 g/mol

TEST	Specification
Appearance (Color)	Yellow to Orange
Appearance (Form) Powder or Crystals	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Nitrogen	9.2- 9.5 %
ICP Major Analysis Confirms Platinum Component	Confirmed
Purity > = 99.9% Based on Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤1,000.0 ppm

Remarks:

Specification Date : 12/13/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.