

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

4-Formylphenylboronic acid -  $\geq 95.0\%$ **Product Number:****431966**

CAS Number:

87199-17-5

MDL:

MFCD00151823

Formula:

C<sub>7</sub>H<sub>7</sub>BO<sub>3</sub>

Formula Weight:

149.94 g/mol

**TEST****Specification**

Appearance (Color)

White to Yellow and Faint Orange to

Orange and Faint Beige to Beige

Appearance (Form)

Conforms to Requirements

Crystals to Crystalline Powder to Powder and/or Chunks

Proton NMR Spectrum

Conforms to Structure

Purity (Titration by NaOH)

 $\geq 95.0\%$ 

Specification: PRD.1.ZQ5.10000021762

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.