


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Silver acetate - ReagentPlus®, 99%

Product Number: 216674
CAS Number: 563-63-3
MDL: MFCD00012446
Formula: C₂H₃AgO₂
Formula Weight: 166.91 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

White to Off-White or Grey

Appearance (Form)

Powder or Chunks

X-Ray Diffraction

Conforms to Structure

Titration by KSCN (% Ag)

63.7 - 65.5 %

Specification: PRD.0.ZQ5.10000022974

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.