


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
trans,trans-Dibenzylideneacetone - 98%

Product Number: 246425
CAS Number: 35225-79-7
MDL: MFCD00004790
Formula: C₁₇H₁₄O
Formula Weight: 234.29 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Yellow to Yellow-Green
Appearance (Form)	Powder or Crystalline Powder or Crystals or Granules
Proton NMR spectra	Conforms to Structure
Carbon	84.9 - 89.3 %
Purity (HPLC)	≥ 97.5 %

Specification: PRD.3.ZQ5.10000043708

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.