

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

EX-CELL® Sp2/0 Serum-Free Medium for Sp2/0 Cells – Animal-component free, Protein free, Chemically defined, without L-glutamine, liquid, sterile-filtered, suitable for cell culture

Product Number: **14660C**

Storage Temperature: 2 - 8 °C

### TEST

### Specification

Appearance (Turbidity)	Clear
Appearance (Color)	Report Only (Color)
Appearance (Form)	Solution
free of ppt	
pH	6.9 - 7.3
Osmolality	305 - 335 mOs/kg
Sterility by USP guidelines	No Microbial Growth Detected
Endotoxin Level	
Record	
Cell Culture Test	≥ 2000000
Growth Promotion/Cytotoxicity	
% Viability	≥ 95 %
Cell Line	mSx15 Cell Line

Specification: PRD.0.ZQ5.10000037489

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.