

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

5-Bromo-4-chloro-3-indolyl β-D-glucuronide cyclohexylammonium salt - tablet

Product Number:

B8049

CAS Number:

114162-64-0

MDL:

MFCD00058543

Formula:


C<sub>14</sub>H<sub>13</sub>BrClNO<sub>7</sub> · C<sub>6</sub>H<sub>13</sub>N

Formula Weight:

521.79 g/mol

Storage Temperature:

-20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Round Tablet
Weight mg	36 - 44
Dissolution Time * One tablet per 20 mL purified water	≤ 25.0 min
Suitability	Suitable
Suitable for detection of gus gene expression; compares to untableted substrate	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.