

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Sodium Channel, Pan antibody produced in mouse - ~1 mg/mL, clone K58/35, purified immunoglobulin

Product Number:

S8809

MDL:

MFCD02097413

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Solution
Protein Concentration by E280	0.9 - 1.3 mg/ml
Immunoblotting:	Conforms
Detects Rat Brain Extract	

Specification: PRD.0.ZQ5.10000028205

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.