


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Bromopyridine-2-carbonitrile - 97%

Product Number: **703966**
CAS Number: 55758-02-6
Formula: C₆H₃BrN₂
Formula Weight: 183.01 g/mol

**TEST****Specification**

Appearance (Color)	White to Tan
Appearance (Form)	Conforms to Requirements
Solid, Powder, Crystals, and/or Chunks	
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃	96.5 - 103.5 %
After Oxygen Combustion	
Purity (GC)	≥ 96.5 %
Recommended Retest Period	-----
24 Months	

Specification: PRD.0.ZQ5.10000031562

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.