

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: **P5976**

Storage Temperature: -20 °C

TEST	Specification
Formulation	Conforms
Lyophilized from 0.2 um-filtered solution in PBS, pH 7.4, with 5% trehalose.	
Bioactivity	Conforms
ND50: 0.1 to 0.5 ug/mL (Cell neutralization, NR6R-3T3 fibroblasts)	
Working Dilution	Conforms
By Indirect Immunoblotting: Working antibody concentration of 0.1 ug/ml detects the recombinant human PDGF-BB protein	
Endotoxin Assay	Conforms
< 0.1 EU/ug	

Specification: PRD.0.ZQ5.10000086332

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.