

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
20S Proteasome Fraction from rabbit - ≥95% (SDS-PAGE), solution

Product Number: **P3988**

MDL: MFCD01091499

Storage Temperature: -70 °C

TEST	Specification
------	---------------

Formulation	Conforms
Solution in 50 mM HEPES, 150 mM NaCl, 1 mM DTT, pH 7.6.	
mg protein/ml	≥ 0.4
Purity by SDS Electrophoresis	≥ 95 %

Specification: PRD.0.ZQ5.10000046080

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.