


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Isopropylaniline - 97%

Product Number: **268887**
CAS Number: 643-28-7
MDL: MFCD00007720
Formula: C₉H₁₃N
Formula Weight: 135.21 g/mol



TEST

Specification

Appearance (Color)
Brown and Colorless to Very Dark Orange
Appearance (Form)
Proton NMR Spectrum
Purity (GC)

Colorless to Very Dark Brown-Yellow and Colorless to Very Dark
Liquid
Conforms to Structure
≥ 96.5 %

Specification: PRD.0.ZQ5.10000029582

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.