

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Nonafluorobutane-1-sulfonic acid - 97%

Product Number: 562629 CAS Number: 375-73-5

MDL: MFCD01320794
Formula: C4HF9O3S
Formula Weight: 300.1 g/mol
Storage Temperature: 2 - 8 ℃

TEST

Appearance (Color)

Appearance (Form)
Infrared spectrum
Fluorine NMR Spectra
Purity (Titration by NAOH)

Specification: PRD.0.ZQ5.10000025060

Specification

Colorless

Liquid Conforms to Structure Conforms to Structure 96.5 - 103.5 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.