


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Clofilium tosylate - >97%, solid

Product Number: **C2365**
CAS Number: 92953-10-1
MDL: MFCD00069233
Formula: C₂₁H₃₇ClN · C₇H₇O₃S
Formula Weight: 510.17 g/mol



TEST

Specification

Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder
Proton NMR spectrum	Conforms to Structure
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
5 mg + 1 mL, MeOH	
Solubility (Color)	Colorless
Solubility (Turbidity)	Slightly Hazy
3 mg + 1 mL, H ₂ O	

Specification: PRD.0.ZQ5.10000058372

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.